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1 General

- The product may be used by trained personnel as well as by lay persons.
- The product is intended to be used to supply medical oxygen for inhalation within homecare oxygen therapy as well as in medical settings.
The instructions of the physician are to be observed.
- The following instructions must be followed by untrained users.
- These instructions reflect the state of the art at the time of publication.
- This product complies with the provisions of the EU Directive 93/42/EEC (Medical Device Directive) and Medical Device Regulation (EU) 2017/745.
- The instructions for use comply with DIN EN 82079-1.

2 Safety and Responsibilities

IMPORTANT! READ CAREFULLY BEFORE USE
KEEP FOR FUTURE REFERENCE

- Please observe the following instructions and care instructions before use.
- If in doubt about the use and connection of the product, please contact customer service.
- If additional oxygen is administered care should be taken to ensure that smoking and other naked flames are forbidden in the vicinity of oxygen-carrying equipment and fittings.
- Do not smoke!
- The product draws air from the immediate surrounding and must therefore not be used in a toxic atmosphere.
- Do not connect the nebuliser system to an endotracheal tube or similar device without an alternate exhalation route.
- Keep the product free of oil and grease, otherwise there is a risk of explosion. This includes the exposure to cleaning alcohols, hand cream and adhesives – in contact with highly compressed oxygen explosive reactions may occur.
- Do not use the product connected to a ventilator or in an MR area.
- Use only oxygen and air.
- A maximum length of 18 m (tubing and inhalation unit) should not be exceeded. Otherwise, the pressure and the oxygen flow cannot be guaranteed.

- The product must not be kinked, stretched or bent.
- Twisting of the tubing should be avoided. Use a swivel connector to prevent twisting.
- **Single use products may not be reused!**
They must be disposed of immediately after use in accordance with chapter 9.
The reuse of single use products represents a potential risk for the patient or the user. It can result in contamination and/or reduced functionality.
- Contamination and/or limited functionality of the product can lead to injury, illness or death of the patient.


3 Modification of the product

- Apart from the described maintenance instructions, the product may not be modified or repaired by the operator!
- The manufacturer's responsibility is limited to the dimensional and physiological suitability of the product.
- Improper use, modifications to the product or negligence will void any warranty. This also applies to unauthorised repairs. Normal signs of wear and tear on the device are also excluded.

4 Before using the product

- Do not expose the product to direct sunlight or other light sources.
Store and transport the product in a dry and dust-free place in its original packaging at -10°C to $+50^{\circ}\text{C}$.
- Remove all delivered parts from the packaging.
If necessary, remove the respective plastic covers and dispose of the packaging according to the current national regulations.
- Before each use verify that parts or components do not show any signs of hardening, damage, breakage or cracks. The product must also not be contaminated or soiled. If any of the conditions described above are present, the device must not be put into operation and must be disposed of in accordance with chapter 9.
- In case of doubt perform a comparison check with a brand-new product.
→ Damaged products must not be used!
After consulting customer service, the product must be disposed of according to chapter 9.
- The product contains ABS, PS, PVC.
- Please note that the product is supplied in a non-sterile condition.

- The hygienic preparation of the product lies within the clinical responsibility of the user! The plastic product can be disinfected with a suitable cleaning agent (e.g. alcohol) before use.
Sterilization is not permitted!
- To attach and detach the tubing, the tubing must be held at the connector to prevent cracks from forming between the connector and the tubing.
- Unscrew the lid of the nebuliser and fill the container with the prescribed medication. Screw the lid back on.
- Attach the mask/mouthpiece securely to the corresponding opening on the container lid.

- Connect one end of the tubing to the nebuliser socket and the other end to the gas source.
- Set the gas supply to the desired flow rate.
Place the mask over mouth and nose onto the face. Pull the elastic over the head and under the ears at the nape of the neck. For a tight fit of the mask, tighten the elastics and adjust the metal nose piece to the patient's face.
- Ensure that aerosol nebulisation is working.
In case of problems with nebulisation of the medication, lightly tap against the nebuliser.
- For total lengths over 15 m, a water trap should be used to collect any condensation.
- When completely assembled, check that all parts are connected tightly.
- The product is now ready for use.

5 Using the product

- The product must not come into contact with sharp objects or rub against an edge during use.
- After successful assembly, the product can be attached to the oxygen source. Please note the requirements of any other connection devices and products.
- Open the gas source and inhale according to the instructions of the physician.
- The flow rate can be set to max. 6 lpm.
- Information on particle size and TOR (Total Output Rate) can be obtained from customer service.

6 Product maintenance / Reprocessing

- Single use products may not be reused!
- The duration of use is limited to a maximum of 7 days.
In case of contamination or defects, replace the product immediately.
- The product must not be used on more than one patient.
- Thoroughly rinse the nebuliser, oxygen mask and headpiece under running water after each use. Use distilled water.
- Dry well and screw back together.
- After termination of the application, the product must be disposed of according to chapter 9.

7 Accessories

- The product can be connected to other tubes via a connector.
- The user is responsible for the correct combination with other medical products.

8 Replacement parts and repairs

- In case of functional or material problems, contact customer service immediately.

9 Disposal

- The product must be disposed of after use.
- Please note the expiry date if applicable.
The product must not be used after the expiry date!
- Safe disposal is carried out in accordance with current national regulations.

11 Symbolerklärung / objaśnienie symboli / Explanation of symbols

		
 Deutschland Kennzeichen	Symbol Niemiec	Symbol for Germany
 Polen Kennzeichen	Symbol Polski	Symbol for Poland
 Groß-Britannien Kennzeichen	Symbol Wielkiej Brytanii	Symbol for Great Britain
 Gebrauchsanweisung beachten	Przestrzegać instrukcji użytkowania	Observe the operating instructions
 Warnhinweis	Wskazówka ostrzegawcza	Warning
 Hersteller	Producent	Manufacturer
 EU-Repräsentant	Przedstawiciel UE	EU Representative
 CE-Symbol	Symbol CE	CE symbol
 Nicht verwenden wenn Packung beschädigt ist	Nie użytkować, jeśli opakowanie jest uszkodzone	Do not use if package is damaged
 Zur einmaligen Verwendung / nur zur Verwendung an einem Patienten	Do jednorazowego użytku / tylko do użytku na jednym pacjencie	Single-use only / Only to be used on one patient
 Medizinprodukt	Wyrób medyczny	Medical Device
 Importeur	Importer	Importer